

03[®] Regional Oximetry

Available for Adult, Paediatric, Infant and Neonatal Applications



The 03 Regional Oximetry platform has been expanded to allow monitoring of infant and neonatal patients <10 kg

- > 03 may help clinicians monitor cerebral oxygenation in situations in which peripheral pulse oximetry alone may not be fully indicative of the oxygen in the brain
- > 03 integrates with Masimo SET[®] pulse oximetry on Root[®], providing clinicians with expanded visibility of a patient's oxygenation status
- > 3% ARMS¹ trending accuracy specification on neonatal patients
- > With its reduced size and flexible design, the 03 neonatal sensor easily conforms to and allows for ergonomic application on small foreheads

O3 Display

Δbase

Displays the difference between current rSO₂ and user-defined baseline

AUC

Area Under the Curve index quantifies the depth and duration of patient stay below user-defined rSO₂ low alarm limit

ΔSpO₂

Displays the difference between SpO₂ (from the Radical-7®, if applicable) and rSO₂



rSO₂
Tissue oxygen saturation


ΔcHbi
Displays an index representing the sum of the ΔO₂Hbi and ΔHHbi components of the rSO₂ calculation

ΔHHbi
Displays an index representing the change in the deoxyhaemoglobin component of the rSO₂ calculation

ΔO₂Hbi
Displays an index representing the change in the oxyhaemoglobin component of the rSO₂ calculation

O3 Monitoring

The Root patient monitoring and connectivity hub offers plug-and-play monitoring with Masimo Open Connect® (MOC-9®) modules.*



Apply the appropriate O3 sensors to the forehead:

- > Infant and Neonatal Adhesive Sensor (<10 kg)
- > Paediatric Adhesive Sensor (≥5 kg and <40 kg)
- > Adult Adhesive Sensor (≥40 kg)



Connect the O3 sensors to an O3 MOC-9 module (up to two sensors per module)



Connect the O3 MOC-9 module to one of three MOC-9 ports on Root

O3 MOC-9 Module Specifications

PHYSICAL CHARACTERISTICS	ENVIRONMENTAL
Length (including cable)..... 12.1 ft (3.7 m)	Operational Temperature..... 32 to 104° F (0 to 40° C)
Width..... 1.8 in (4.6 cm)	Storage Temperature..... -40 to 158° F (-40 to 70° C)
Thickness..... 0.6 in (1.5 cm)	Operating and Storage Humidity..... 10 to 95%, non-condensing
Weight..... 7.1 oz max (200 g max)	Altitude..... Up to 12,000 ft (3700 m)

O3 Sensor Specifications

	ENVIRONMENTAL
Application Site..... Forehead	Operating Temperature at Ambient Humidity..... 41 to 104° F (5 to 40° C)
Wavelengths..... 4	Storage Temperature at Ambient Humidity..... -40 to 140° F (-40 to 60° C)
Neonatal rSO₂ Sensor Accuracy (ARMS)² <10kg	Storage Humidity..... 15% to 95%, 86 to 140° F (30 to 60° C)
Trending Regional Oxygen Saturation (rSO ₂)..... 3%	
Paediatric rSO₂ Sensor Accuracy (ARMS)² ≥5 kg and <40 kg	
Absolute Regional Oxygen Saturation (rSO ₂)..... .5%	
Trending Regional Oxygen Saturation (rSO ₂)..... .3%	
Adult rSO₂ Sensor Accuracy (ARMS)² ≥40 kg	
Absolute Regional Oxygen Saturation (rSO ₂)..... .4%	
Trending Regional Oxygen Saturation (rSO ₂)..... .3%	

¹ ARMS accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within ± ARMS of the reference measurements in a controlled study.

* In countries with regulatory approval and Root devices with the correct software version.

The O3 System with infant and neonatal sensors is not licensed for sale in Canada.

For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo U.S.
Tel: 1 877 4 Masimo
info-america@masimo.com

Masimo International
Tel: +41 32 720 1111
info-international@masimo.com

