

# Radius-7™

## The Power of Masimo's Breakthrough Measurements in a Patient-worn Monitor

### > Breakthrough measurements

- Oxygen saturation (SpO<sub>2</sub>) and pulse rate monitoring with Masimo SET® Measure-through Motion and Low Perfusion™ pulse oximetry for reduced false alarms and increased true alarm detection during challenging conditions<sup>1,2</sup>
- Respiration rate monitoring with rainbow Acoustic Monitoring™ for increased patient tolerance compared to capnography<sup>3</sup>



### > Small, lightweight, and wearable for untethered monitoring and ambulation

- Designed to promote greater patient comfort and independence
- Reduces the need for nurses to disconnect the monitor each time the patient gets out of bed

### > Standard wireless communication

- Short-range communication to Root® with secured Bluetooth
- Upgradeable to WiFi for direct communication throughout the hospital to the Patient SafetyNet™\*\* remote monitoring system



# FLEXIBLE FUNCTIONALITY

- > Parameter and waveform views facilitate on-the-spot assessment of patient status
- > Configurable display shut-off minimizes patient distraction and supports privacy while ambulating
- > Alarms through Root when in range or directly on Radius-7 when out of range
- > Two “hot-swappable,” rechargeable modules (one on Radius-7, one charging in Root) minimize monitoring interruption
- > 12-hour battery life facilitates just one module exchange per shift
- > Single-patient use armband minimizes cross-contamination risk while ensuring patient comfort
- > Seamless module pairing with Root to facilitate nursing workflows



Each Radius-7 comes with two rechargeable, “hot-swappable” modules.

## SPECIFICATIONS

### PHYSICAL CHARACTERISTICS

Dimensions ..... 5.1" x 2.8" x 1.2" (130 mm x 70 mm x 30 mm)  
 Weight ..... 0.34lbs. (155g)

### INDICATORS

SpO<sub>2</sub>, pulse rate, Respiratory Rate (RRa), Pleth Variability Index (PVI®), perfusion index, pleth waveform, RRa waveform, alarm status  
 Display Update Rate ..... 1 second

### ALARMS

Audible alarms sound on Root when connected via Bluetooth. Audible alarms sounds on Radius-7 when not connected to Root via Bluetooth. Visual alarms are displayed on both Radius-7 and Root.

### EACH RADIUS-7 INCLUDES

- 2 Battery modules
- 1 Instrument module
- 1 Battery charging adapter

### BATTERY MODULE

Type ..... Lithium ion  
 Capacity (battery life) ..... 12 hours  
 Charging Time ..... ≤ 6 hours

### COMMUNICATION

Radio ..... Bluetooth  
 WiFi ..... 802.11 A/B/G<sup>4</sup>

### Communication Compliance

USA ..... FCC ID: VKF-RADIUS  
 Model: RADIUS  
 FCC parts 15.247  
 Canada ..... IC: 7362A-RADIUS  
 IC Model: RADIUS  
 RSS-210  
 Europe R&TTE Directive ..... EN 300 328, EN 301 489-17

For complete specifications, including measurements, see Operator’s Manual.

<sup>1</sup> Taenzer A.H. et al. *Anesthesiology*. 2010 Feb;112(2):282-7. <sup>2</sup> Shah N. et al. *J Clin Anesth*. 2012 Aug;24(5):385-91. <sup>3</sup> Patino M et al. *Paediatr Anaesth*. 2013 Dec;23(12):1166-73. <sup>4</sup> Not available in all models.

\* The use of the trademarks Patient SafetyNet and PSN is under license from University HealthSystem Consortium.

Caution: Federal law restricts this device to sale by or on the order of a physician.

**Masimo U.S.**  
 Tel: 1 877 4 Masimo  
 info-america@masimo.com

**Masimo International**  
 Tel: +41 32 720 1111  
 info-international@masimo.com

