

| Newborn Sensors

Featuring Masimo SET® Measure-through Motion and Low Perfusion™ Pulse Oximetry




Newborn Sensors automatically configure Masimo SET® and rainbow SET™ devices for the fastest averaging time and maximum sensitivity settings




Newborn Sensor Options




- > VelAid™ SofTouch design allows for quick, easy sensor application to and repositioning on newborn skin
- > Hook and loop attachment strip keeps the sensor secure, even when site is wet



RD SET™ SENSORS



LNCS® SENSORS



M-LNCS™ SENSORS

Specifications

ACCURACY (ARMS) ¹	WEIGHT RANGE
Oxygen Saturation (%SpO ₂) Accuracy Range 70-100%	RD SET Newborn Neonatal <3 kg, hand or foot application
No Motion Pediatrics/Infants 2%	RD SET Newborn Infant/Pediatric 3-10 kg, thumb or great toe application
No Motion Neonates 3%	10-30 kg, finger or toe application
Motion Pediatrics/Infants/Neonates 3%	LNCS Newborn Neonatal <3 kg, hand or foot application
Low Perfusion Pediatrics/Infants 2%	LNCS Newborn Infant/Pediatric 3-10 kg, thumb or great toe application
Low Perfusion Neonates 3%	10-30 kg, finger or toe application
Pulse Rate (bpm) Accuracy Range 25-240 bpm	M-LNCS Newborn Neonatal <3 kg, hand or foot application
No Motion Pediatrics/Infants/Neonates 3 bpm	M-LNCS Newborn Infant/Pediatric 3-10 kg, thumb or great toe application
Motion Pediatrics/Infants/Neonates 5 bpm	10-30 kg, finger or toe application
Low Perfusion Pediatrics/Infants/Neonates 3 bpm	
INFORMATION	ORDERING INFORMATION
Single-patient-use / Non-sterile / Packaged 20 per box / Does not contain natural rubber latex	RD SET Newborn Neonatal 4013
	RD SET Newborn Infant/Pediatric 4012
	LNCS Newborn Neonatal 2412
	LNCS Newborn Infant/Pediatric 2413
	M-LNCS Newborn Neonatal 2519
	M-LNCS Newborn Infant/Pediatric 2520

¹ ARMS Accuracy is a statistical calculation of the difference between device measurements and reference measurements. Approximately two-thirds of the device measurements fell within ± ARMS of the reference measurements in a controlled study.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo U.S.
Tel: 1 877 4 Masimo
info-america@masimo.com

Masimo International
Tel: +41 32 720 1111
info-international@masimo.com

