







Sensor & Cable Compatibility

Technology ▶	Masimo SET® or Masimo rainbow SET™		
	Cable ▶	M-LNC/rainbow® RC	RD SET™/RD rainbow SET™
Connector ▶	DB9	M15	RD15
			
Sensor			
LNCS® Disposables	Y	N	N
LNCS Reusables	Y	N	N
M-LNCS™ Disposables	N	Y	N
M-LNCS Reusables	N	Y	N
RD SET Disposables	N	N	Y
RD SET Reusables	N	N	Y

Minimum Software Compatibility Requirements*

Parameters ▶	SET® Parameters					
	Cable ▶	LNC	Red LNC	M-LNC	rainbow® RC	RD SET
Connector ▶	DB9		M15		RD15	
Sensor						
LNCS Disposables	Compatible	7.0	Not compatible	Not compatible	Not compatible	Not compatible
LNCS Reusables	Compatible	7.0	Not compatible	Not compatible	Not compatible	Not compatible
M-LNCS Disposables	Not compatible	Not compatible	Compatible	7.1	Not compatible	Not compatible
M-LNCS Reusables	Not compatible	Not compatible	Compatible	7.1	Not compatible	Not compatible
RD SET Disposables	Not compatible	Not compatible	Not compatible	Not compatible	Compatible	7.1
RD SET Reusables	Not compatible	Not compatible	Not compatible	Not compatible	Compatible	7.1
Trauma and Newborn	4.1	7.0	4.1	7.1	4.1	7.1

Compatible
 Not compatible

* Masimo technology board software, sensor, and sensor cable combinations have been validated and tested for compliance per ISO 80601-2-61.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo Americas
 tel 1-877-452-7455
 info-america@masimo.com

Masimo International
 tel +41-32-720-1111
 info-international@masimo.com

