

rainbow® Sensor & Cable Compatibility

		rainbow SET™		SET®		
		RD rainbow SET™	rainbow® RC	RD SET™	M-LNC	Red LNC/LNC
Masimo Technology ▶						
rainbow® Sensor ▼	Patient Cable ▶					
	Connector ▶					
		RD15	M15	RD15	M15	DB9
RD rainbow SET-2 Disposables		Y	N	Y	N	N
rainbow® Disposables R		N	Y	N	Y	N
rainbow® Disposables R1		N	Y	N	Y	N
rainbow ReSposable® R2		N	Y	N	Y	N
rainbow® SpHb® Spot-check Reusables		N	Y	N	N	N
rainbow® SpCO® Reusables		N	Y	N	Y	N

Minimum Software Compatibility Requirement*

		RD rainbow SET™					RD SET™	
Patient Cable ▶	rainbow® Sensor ▼	SpHb®	SpMet®	SpCO®	PVi®	SET®	PVi	SET®
	RD rainbow SET-2 Disposables	7.4	7.1		4.8	No Min	4.8	No Min

		rainbow® RC					M-LNC	
Patient Cable ▶	rainbow® Sensor ▼	SpHb	SpMet	SpCO	PVi	SET®	PVi	SET®
	rainbow® Disposables R		7.1	7.1	4.8	No Min	4.8	No Min
	rainbow® Disposables R1	7.4	7.1		4.8	No Min	4.8	No Min
	rainbow ReSposable® R2	7.4	7.1		4.8	No Min	4.8	No Min
	rainbow® SpHb Spot-check Reusables**	7.4	7.4		7.4	7.4		
	rainbow® SpCO Reusables		7.1	7.1	4.8	No Min	4.8	No Min

* Masimo technology board software, sensor, and sensor cable combinations have been validated and tested for compliance per ISO 80601-2-61.

** Device must have SpHb parameter enabled.

RD rainbow SET Sensors are not licensed for sale in Canada.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo Americas
tel 1-877-452-7455
info-america@masimo.com

Masimo International
tel +41-32-720-1111
info-international@masimo.com

