

SET® and rainbow® Sensor & Cable Compatibility

Masimo Technology ▶		rainbow SET®		SET®			
Patient Cable ▶		RD rainbow SET®	rainbow® RC	RD SET®	M-LNC	Red LNC/LNC	
Sensor ▼		Connector ▶	RD15	M15	RD15	M15	DB9
RD rainbow SET-2 Disposables			Y	N	Y	N	N
RD rainbow 8λ SpCO® Disposables			Y	N	Y	N	N
RD rainbow™ 4λ Disposables			Y	N	Y	N	N
LNCS-II™ rainbow 8λ SpHb® Spot-check Reusables			Y	N	Y	N	N
LNCS-II rainbow 8λ SpCO Reusables			Y	N	Y	N	N
RD SET Disposables & Reusables			Y	N	Y	N	N
Sensor ▼							
rainbow® Disposables R			N	Y	N	Y	N
rainbow® Disposables R1			N	Y	N	Y	N
rainbow® SpHb Spot-check Reusables			N	Y	N	N	N
M-LNCS® Disposables & Reusables			N	Y	N	Y	N
Sensor ▼							
LNCS® Disposables & Reusables			N	N	N	N	Y

RD rainbow 4λ and RD rainbow 8λ SpCO sensors are not licensed for sale in Canada.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo Americas
tel 1-877-452-7455
info-america@masimo.com

Masimo International
tel +41-32-720-1111
info-international@masimo.com



Minimum Software Compatibility Requirement*

Patient Cable ▶		RD rainbow SET®					RD SET†			
Sensor ▼	Parameter ▶	SpHb [®] /SpOC™	SpMet®	SpCO®	PVi®	RRp®	SET®	PVi	RRp†	SET®
RD rainbow SET-2 Disposables		7.4	7.1		7.1	7.9	7.1	7.1	7.9	7.1
RD rainbow 8λ SpCO Disposables			7.1	7.1	7.1	7.9	7.1	7.1	7.9	7.1
RD rainbow 4λ Disposables					7.1	7.9	7.1	7.1	7.9	7.1
LNCS-II rainbow 8λ SpHb Reusables		7.4	7.4		7.4	7.9	7.4			
LNCS-II rainbow 8λ SpCO Reusables			7.1	7.1	7.1	7.9	7.1	7.1	7.9	7.1
RD SET Disposables & Reusables					7.1	7.9	7.1	7.1	7.9	7.1

Patient Cable ▶		rainbow® RC					M-LNC			
Sensor ▼	Parameter ▶	SpHb/SpOC	SpMet	SpCO	PVi	RRp	SET®	PVi	RRp	SET®
rainbow® Disposables R1		7.4	7.1		7.1	7.9	7.1	7.1	7.9	7.1
rainbow® Disposables R			7.1	7.1	7.1	7.9	7.1	7.1	7.9	7.1
rainbow® SpHb Spot-Check Reusables††		7.4	7.4	--	7.4	7.9	7.4			
rainbow® SpCO Reusables			7.1	7.1	7.1	7.9	7.1	7.1	7.9	7.1
M-LNCS Disposables & Reusables					7.1	7.9	7.1	7.1	7.9	7.1

Patient Cable ▶		Red LNC					LNC			
Sensor ▼	Parameter ▶	SpHb/SpOC	SpMet	SpCO	PVi	RRp	SET®	PVi	RRp	SET®
LNCS Disposables & Reusables					7.1	7.9	7.1	7.1	7.9	7.1

* Masimo technology board software, sensor, and sensor cable combinations have been validated and tested for compliance per ISO 80601-2-61. Minimum software versions for MX technology boards shown. Additional compatibility for MS technology boards for all sensor and cable combinations listed is as follows: SET - no minimum, PVi 4.8.

† RRp is automatically provided on all versions of Rad-G.

†† Device must have SpHb parameter enabled.