




NomoLine® Device Compatibility Chart

The Masimo NomoLine family of sampling lines is compatible with many existing capnography and gas monitoring devices.

Masimo Sampling Line	Compatible Devices with Examples
 <p>NomoLine</p>	<p>Masimo devices such as:</p> <ul style="list-style-type: none"> > Root® with ISA™ > Rad-97™ with NomoLine Capnography
 <p>NomoLine-O™</p>	<p>Medtronic (Oridion) Microstream Devices</p> <ul style="list-style-type: none"> > Capnostream 35 > Microstream EtCO2 Modules
 <p>NomoLine with Luer</p>	<p>Devices with Water Traps</p> <ul style="list-style-type: none"> > GE with Datex Ohmeda > Dräger

NOTE: Devices without NomoLine compatibility include Philips Respironics, Nihon Kohden, and Nonin

The GE®, MINDRAY®, DRÄGER®, PHILIPS®, Philips Respironics®, NIHON KOHDEN®, and NONIN® trademarks are property of their respective owners and are used for identification purposes only. Use of these names, trademarks, or brands does not imply any affiliation with or endorsement.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions. For professional use. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo Americas
tel 1-877-452-7455
info-america@masimo.com

Masimo International
tel +41-32-720-1111
info-international@masimo.com

