

Bridge™

Drug-Free Opioid Withdrawal Device



Bridge is the first FDA-cleared, drug-free, non-surgical device to use neuromodulation to aid in the reduction of symptoms associated with opioid withdrawal.

In a study of 73 adult patients with opioid use disorder, a reduction in opioid withdrawal symptoms was observed when Bridge was used during active withdrawal. The study demonstrated Bridge to be:

FAST

- > Provides symptom relief in as little as 20 minutes¹


EFFECTIVE

- > May reduce opioid withdrawal symptoms
- > Helps users successfully make treatment progress:
 - 88% of patients successfully transitioned into a medically assisted therapy (MAT) program¹

SAFE²

- > No reported side effects with minimal risks during device placement
 - < 0.17% Incidence of significant pain at the lead implantation site
 - < 0.91% Incidence of bleeding at the lead implantation site
 - < 0.91% Incidence of dermatitis at the lead implantation site
 - 0% Incidence of syncope

Mean Clinical Opioid Withdrawal Scale (COWS) Scores While Using Masimo Bridge During Withdrawal



What is Bridge?


Bridge is a small electrical nerve stimulator device that contains a battery-powered chip and wires that are applied around a patient's ear. It requires a prescription and is offered to qualified healthcare professionals with training. The Bridge kit includes everything needed to apply the device.

Applied by a qualified healthcare professional in a non-surgical, in-office procedure

Fits comfortably behind your ear

Each device provides five days of continuous relief

Sends gentle electrical impulses through wires to nerves around your ear



Stimulated nerves transmit these impulses to the brain, reducing withdrawal symptoms

Who Can Bridge Help?

Bridge can be used for patients* experiencing opioid withdrawal symptoms, while undergoing treatment for opioid use disorder, when:

- Initiating treatment
- Transitioning to naltrexone
- Tapering off medication-assisted treatment

Specifications

Ordering Information

Masimo Part Number 4790

* Contraindications: Use of cardiac pacemakers, hemophilia, psoriasis vulgaris.

¹ Miranda A et al. *Am J Drug Alcohol Abuse*. 2018;44(4):498]. *Am J Drug Alcohol Abuse*. 2018;44(1):56-63. doi:10.1080/00952990.2017.1295459. ² Roberts A et al. Minimal adverse effects profile following implantation of periauricular percutaneous electrical nerve field stimulators: a retrospective cohort study. *Medical Devices (Auckland, N.Z.)* vol. 9 389-393. 3 Nov. 2016, doi:10.2147/MDER.S107426.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, and precautions.

Masimo U.S.
Tel: 1 877 4 Masimo
info-america@masimo.com

Masimo International
Tel: +41 32 720 1111
info-international@masimo.com

